

# **CODING OF A SYSTEM**

## **I. MICROCONTROLLER PROGRAMMING (ARDUINO IDE)**

```
#include <ESP8266WiFi.h>
char ssid[] = "Project";
char pass[] = "12345678";
WiFiClient client;
#include "ThingSpeak.h"
//#include <WiFi.h>
//#include "WiFi.h"
#define SECRET_CH_ID 1736899
#define SECRET_WRITE_APIKEY "OEMFMV390IANB4YF"
#include <LCD_I2C.h>
float x=0.0;
int y=0;
LCD_I2C lcd(0x27);
#include <Wire.h>
#include "MAX30100_PulseOximeter.h"
float vref = 3.3;
float resolution = vref/1023;
#define REPORTING_PERIOD_MS    1000
PulseOximeter pox;
int keyIndex = 0;
String myStatus = "";
unsigned long myChannelNumber = SECRET_CH_ID;
const char * myWriteAPIKey = SECRET_WRITE_APIKEY;
uint32_t tsLastReport = 0;

void onBeatDetected()
{
```

```

    Serial.println("Beat!");
}

void setup()
{
    Serial.begin(115200);
    lcd.begin();
    lcd.begin(false)

    lcd.backlight();
    lcd.setCursor(0, 0);
    lcd.print("PATIENT");
    lcd.setCursor(0, 1);
    lcd.print("MONITORING");
    delay(5000);
    lcd.clear();
    if(WiFi.status() != WL_CONNECTED){

    while(WiFi.status() != WL_CONNECTED){
        WiFi.begin(ssid, pass); // Connect to WPA/WPA2 network. Change this
line if using open or WEP network
        // Serial.print(".");
        delay(5000);
        break;
    }
    // Serial.println("\nConnected.");
}

    WiFi.mode(WIFI_STA);
    ThingSpeak.begin(client);
    Serial.print("Initializing pulse oximeter..");

```

```

    if (!pox.begin()) {
        Serial.println("FAILED");
        for(;;);
    } else {
        Serial.println("SUCCESS");
    }
    pox.setOnBeatDetectedCallback(onBeatDetected);

}

void loop()
{
    pox.update();
    if (millis() - tsLastReport > REPORTING_PERIOD_MS) {
        x=pox.getHeartRate();
        y=pox.getSpO2();
        Serial.print("Heart rate:");
        Serial.print(x);
        Serial.print("bpm / SpO2:");
        Serial.print(y);
        Serial.println("%");

        if((x>70)&&(y>90))
        {

            lcd.setCursor(0, 0);
            lcd.print("Heart:");
            lcd.print(x);
            lcd.setCursor(0, 1);
            lcd.print("SP02:");

```

```
    lcd.print(y);
    delay(2000);
    lcd.clear();
    float temperature = analogRead(A0);
    temperature = (temperature*resolution);
    temperature = temperature*100;
    Serial.println(temperature);
    lcd.setCursor(0, 0);
        lcd.print("TEM:");
        lcd.print(temperature);
        delay(2000);
        lcd.clear();
        ThingSpeak.setField(1, x);
        ThingSpeak.setField(2,y);
        ThingSpeak.setField(3,temperature);
        int x = ThingSpeak.writeFields(myChannelNumber, myWriteAPIKey);
        delay(60000);
    }
    tsLastReport = millis();

    }

}
```

## II. DATASET TRAINING AND TESTING CODE:

```
import numpy as np
import pandas as pd
from sklearn.model_selection import train_test_split
import matplotlib.pyplot as plt
import seaborn as sns
import pickle

data = pd.read_csv('/content/data.csv')
data.head()
data.shape
X = data.iloc[:, :-1]
X.head()
y = data.iloc[:, -1]
y.head()
data['target'].value_counts()
X_train, X_test, y_train, y_test =
train_test_split(X, y, test_size=0.2, random_state=1)
sns.countplot(x='target', data=data)
plt.show()
X_train.shape
X_train.head()
y_test.shape
y_test.head()

from sklearn.neighbors import KNeighborsClassifier

model = KNeighborsClassifier(n_neighbors=3)
model.fit(X_train, y_train)
filename = 'Patient_model.sav'
```

```

pickle.dump(model, open(filename, 'wb'))
y_pred = model.predict(X_test)
from sklearn import metrics
acc=(metrics.accuracy_score(y_pred,y_test))
print("Accuracy is:",acc)
print("Confusion Matrix is: ",metrics.confusion_matrix(y_pred,y_test))

```

### ***III.DATA CLASSIFICATION AND PREDICTION***

```

import pickle
import urllib.request
import json
from time import sleep

conn =
urllib.request.urlopen("https://api.thingspeak.com/channels/1736899/feeds.json?
results=1")

response = conn.read()
print ("http status code=%s" % (conn.getcode()))
data=json.loads(response)
x=int(data['feeds'][0]['entry_id'])
y=x
conn.close()

while x==y:
    conn =
urllib.request.urlopen("https://api.thingspeak.com/channels/1736899/feeds.json?
results=1")

    response = conn.read()
    #print ("http status code=%s" % (conn.getcode()))
    data=json.loads(response)
    y=int(data['feeds'][0]['entry_id'])
    conn.close()

```

```

conn =
urllib.request.urlopen("https://api.thingspeak.com/channels/1736899/feeds.json?
results=1")
response = conn.read()
print ("http status code=%s" % (conn.getcode()))
data=json.loads(response)
a=float(data['feeds'][0]['field1'])
b=float(data['feeds'][0]['field2'])
c=float(data['feeds'][0]['field3'])

conn.close()
filename = 'Patient_model.sav'
loaded_model = pickle.load(open(filename, 'rb'))
person_reports = [[a,b,c]]
disease_predicted = loaded_model.predict(person_reports)
print("ANALYSING....")
sleep(15)
disease_predicted[0]=0
if disease_predicted[0]==1:
    print("The person may have no disease")
    #sleep(30)
    conn =
urllib.request.urlopen("https://api.thingspeak.com/update?api_key=OEMFMV3
90IANB4YF&field7=NO_DISEASES")
elif disease_predicted[0]==0:
    print("The person may be in Fever take Paracetamol")
    #sleep(30)
    conn =
urllib.request.urlopen("https://api.thingspeak.com/update?api_key=OEMFMV3
90IANB4YF&field7=PARACETAMOL")
elif disease_predicted[0]==2:

```

```
print("The person may be in Hypertension take nisoldipine")
#sleep(30)

conn =
urllib.request.urlopen("https://api.thingspeak.com/update?api_key=OEMFMV3
90IANB4YF&field7=NISOLDIPINE")

elif disease_predicted[0]==3:
    print("The person may have Covid Visit Hospital" )

    conn =
urllib.request.urlopen("https://api.thingspeak.com/update?api_key=OEMFMV3
90IANB4YF&field7=VISIT_HOSPITAL")
```